

Componentes del kit
Components of the kit

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ENGLISH

INSTRUCTIONS FOR FREEZING
SEMEN AT HOME

WARNING:

The cryoprotectant liquid should be put in the home refrigerator as soon as the KIT is received.

Components of the kit

- A) 1 sterile container for the semen sample.
- B) 1 container with cryoprotectant liquid.
- C) 2 pipettes.
- D) 1 sterile gauze.
- E) 4 cryotubes or vials (small containers for cryopreservation).
- F) 1 permanent marker.
- G) 1 tank for transporting biological samples with controlled temperature and a rod inside. It contains a safety tie in the packaging.

Instructions for obtaining the sample

1. Abstain from ejaculation for 3 to 5 days.
2. Wash hands and genitals well with soap and water and rinse with abundant water.
3. Obtain the semen sample by masturbation.
4. Deposit all the material in the sterile container (A).
5. Fill in the questionnaire.


Instructions for freezing the sample

1. Take the cryoprotectant (B) out of the refrigerator before obtaining the sample.
2. Place all the components of the kit on the sterile gauze (D), where the procedure will be carried out. **It is highly important** that pipettes (C) do not come into contact with any object or surface.
3. Mark the four cryotubes (E) with the name, surname and date the sample was obtained (mark each tube).
4. Add the cryoprotectant to the sample: with a pipette, draw up 2 ml of the cryoprotectant liquid. Drop by drop, very slowly, pour the cryoprotectant into the container with the semen sample. Move in circles to mix it well.
5. Distribute the mixture in the 4 cryotubes with the help of the other pipette. Any surplus material can be disposed of as domestic waste.
6. Leave it to rest for 10 minutes.

Instructions for placing the cryotubes in the tank

1. Open the tank (G): a little vapour may come out, this is not dangerous or corrosive.
2. Place the cryotubes in the rod (which is inside the tank).
3. Insert the rod with the cryotubes into the tank.
4. Close the tank with the cover and the tank packaging with the safety tie.
5. To end the process, call the number specified in the transport instructions for the tank to be collected. The tank must be collected within 48 hours.

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